PREVENTIS *SmarTest*®

Calprotectin Home

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Literatur

1. Kawashima K, Ishihara S, Yuki T et al. Fecal calprotectin level correlated with both endoscopic severity and disease extent in ulcerative colitis. *BMC Gastroenterology* (2016) 16:47
2. Najari A, Hiraoka S, Takahashi S. et al. Simultaneous Measurements of Faecal Calprotectin and the Faecal Immunochemical Test in Quiescent Ulcerative Colitis Patients Can Stratify Risk of Relapse. *Journal of Crohn's and Colitis* (2018), 71–76
3. Amcoff K, Stridsberg M, Lampinen M et al. Clinical implications of assay specific differences in f-calprotectin when monitoring inflammatory bowel disease activity over time. *Scandinavian Journal of Gastroenterology* (2017), 52:3, 344-350
4. Matsura R., Watanabe O, Nakamura M. et al. Fecal calprotectin reflects endoscopic activity in patients with small-bowel Crohn's disease according to double-balloon endoscopy findings. *Nagoya J. Med. Sci.* 80 (2018). 257–266

PREVENTIS SmarTest®

Calprotectin Home

REF PST11001EP01

UDI-DI 04050598007697

PREVENTIS

Intended purpose

SmarTest Calprotectin Home is a rapid test, which is performed manually. It includes a smartphone app developed for the quantitative detection of inflammation marker, Calprotectin, in human stool. This test is an in vitro diagnostic medical device for single use only and suitable for self-testing by non-professionals from 18 years old. It serves as an aid in follow up and therapy monitoring of inflammatory bowel disease (e.g. Crohn's disease or ulcerative colitis).

The test is single-use-only.

General information

- **SmarTest Calprotectin Home** is intended for in-vitro diagnostics only. For external use only.
- One test kit contains materials for the performance of **1 test**.
- The performance of the test does not require any special expertise. It can be carried out at home and requires at least 20 minutes.
- Before carrying out the test, carefully read the instructions for use. Follow the instructions for use precisely to set up the application, extract and measure the sample.

- Ensure that your smartphone is sufficiently charged (required battery level of at least 20 %).
- **SmarTest Calprotectin app** requires an active Internet connection in regular intervals, at least every 48 hours to update the latest calibration data.
- Do not switch off the phone or ring tone/alerts while using the product as the timer has an alarm.
- Do not mix reagents that have different batch numbers.
- The test result can be evaluated **no earlier than 15 minutes** and **no later than 17 minutes** after the sample was placed onto the test cassette. **The test can no longer be evaluated 2 minutes after the evaluation period expires.**
- For privacy reasons, your data is only transmitted after encryption. Please find more details on this topic at **Preventis GmbHs Data Privacy Policy** at privacy.preventis.com/smarest/calprotectin and terms.preventis.com.

Safety notes

- Do not drink the contents of the sample collection tube and avoid direct contact with the liquid. In case of skin contact, immediately wash contaminated areas thoroughly.
- People under 18 years old, elderly and people with disabilities must be assisted with setting up the app, extracting the sample, and implementing and evaluating the test.
- Treat the sample as potentially infectious: do not eat, drink, or smoke while handling the sample. We recommend wearing disposable gloves and to always wash/disinfect your hands after testing.
In case of serious incidents, please contact Preventis GmbH immediately at info@preventis.com.
- **Caution:** Please note that children under the age of 4 often naturally have high fecal calprotectin values.

Storage, stability and disposal

- Store the test at room temperature or in the refrigerator (+4 to 30°C). Do not freeze.
- The test is sensitive to humidity, direct sunlight and high temperatures. For this reason, the test must be protected against heat and the **test cassette must be used immediately after opening the packaging (within 5 minutes)**.
- Do not use the test if the packaging or included materials are damaged.
- Do not use the test after the expiration date.
- After use, dispose of all test components with your household waste, including camera test card and test instructions. The paper stool catcher can be flushed down the toilet.

Medical background

Calprotectin is a calcium-binding protein that is produced and released in large quantities by neutrophil granulocytes (neutrophils) in inflammatory processes. In case of inflammatory bowel disease, neutrophils move into the intestinal lumen, where they release calprotectin, which can then be found in the stool. However, based on clinical symptomatology, it is often difficult to distinguish patients with chronic inflammatory bowel diseases (IBD) such as Crohn's disease and ulcerative colitis from those with common irritable bowel syndrome. The additional determination of fecal calprotectin, then, makes it possible to diagnose intestinal inflammations with clinical safety. An increased concentration of fecal calprotectin provides information about the presence of inflammation in the bowels. Non-invasive calprotectin testing enables the physician to spare patients with low concentrations and common inflammatory bowel disease from an invasive colonoscopy, which is

oftentimes unnecessary. Since scientific studies have proven that the fecal calprotectin concentration correlates well with histological and endoscopic results of the disease activity in IBD patients, fecal calprotectin concentration can also be used diagnostically to objectively document therapeutic success. Furthermore, it makes it possible to monitor IBD patients without symptoms to determine and reduce a relapse (recurrence) of IBD early on. In addition, fecal calprotectin serves as a biological marker for inflammatory tumorous changes in the gastrointestinal tract (see **References**, page 88).

Test principle

SmarTest Calprotectin Home is an antibody based strip test for the detection of calprotectin in a stool samples (15 mg). The stool sample is added and dissolved in an extraction solution using a specific sample extraction system. Afterwards, **3 drops** of the solution are placed on the **sample application window** of the rapid test. Existing calprotectin reacts in the test strip with calprotectin antibodies, which are marked with gold. If carried out correctly, a red test line becomes visible after no longer than 15 minutes. The appearance of the control line (C) indicated whether the test functioned correctly (see **Evaluation**, page 86). Using the **SmarTest Calprotectin app**, the quantitative result is then evaluated based on the intensity of the color in the control and test lines. Step-by-step descriptions of the extraction and processing of the sample before measuring, as well as the performance and evaluation of the test are described in the chapters **Sample Extraction** (see page 78), **Test Procedure** (see page 82) and **Evaluation** (see page 86).

Included materials

1 test cassette **TEST**, individually wrapped with desiccant bag *

1 sample collection tube with buffer solution


TUBE

1 camera test card **CARD**

1 paper stool catcher

* To prevent the test strip from absorbing moisture, the desiccant bag is included in the packaging. Dispose of desiccant bag after opening.

Additional required materials

Disposable gloves, a smartphone with **SmarTest Calprotectin** app and an internet connection are needed.

Setting up the SmarTest Calprotectin App

Download and Installation

Attention: The successful installation (app symbol appears on your smartphone) of the **SmarTest Calprotectin** app is required for the evaluation of the rapid test!

The SmarTest Calprotectin app is available for:

✓ iOS, available on the App Store

✓ Android, available on the Google Play Store

You can find the smartphone types and operating systems previously tested by us at calprotectin.preventis.com. If your smartphone is not listed, you can qualify it for the test using the camera test card (CARD). **Attention:** Unfortunately, not all smartphones meet the app's requirements. If you have questions, please contact us at info@preventis.com. If you do not receive the request for device qualification, move to the next step, the **camera test**, page 76, as your smartphone type has then already been approved by us.

mySmarTest

This is a portal for the storage of your test results. In addition, you can use mySmarTest to share test results with medical professionals under „**release results**“ (see **Releasing test results** page 72). You can set up your account when you set up the app, or via „**mySmarTest**“ on the main screen.

Account information


The identity consists of an anonymous, randomly generated user name and a password selected by you. This identity is stored **independently of your test results**. If you forget your user name or password, your account can only be restored, if you saved an email address. You can store the email address optionally upon registration, or under settings in „**mySmarTest**.“

Releasing test results

You need a **patient information form** to provide medical professionals with access to your test results. You will receive this form from your attending physician's office. Open the **SmarTest Calprotectin App** and then press „**mySmarTest**“ on the main screen. Afterwards, you can scan the

QR Code on your patient information form via „Release results“ and thereby share your test results with the medical professionals providing your care. You can manage your access permissions at anytime under „Manage releases.“

Device Qualification

Step 1

Turn the camera test card (CARD) in such a way that the blue stripe is on the right side.

Note: If the app does not request the implementation of a device qualification, please move directly to the **Camera Test** chapter on page 76.

Note: If your device qualification fails repeatedly, please check your camera lens for scratches and stains and remove them, if possible. Also check the CARD for potential damages.


Attention: Make sure that the lighting is even. Avoid shadows and direct sunlight.

If the check continues to fail, please contact us at info@preventis.com.

Step 2

Align the stencil on your display with the six quadrants on the CARD. If aligned correctly, the stencil will turn green. Be sure to hold your smartphone parallel to the CARD until you receive a success message.

Camera Test

Step 1

Place the CARD in front of you with the blue stripe facing up and the numbers showing correctly. Align the stencil on your phone display with the illustration on test cassette ① until you receive a success notification.

Step 2

Repeat this step with the illustration on test cassette ②.

Attention: Carry out the camera test in the same location in which you will later conduct the rapid test to account for the local light conditions. Make sure that the lighting is even. Avoid shadows and direct sunlight.

Note: Be sure to hold your smartphone parallel to the CARD. If the camera test fails repeatedly, please check your camera lens for scratches and stains and remove them if possible. Also check the CARD for potential damages.

If the check continues to fail, please contact us at info@preventis.com.

Sample Extraction

Open the SmarTest Calprotectin app on your smartphone and press the button „Prepare Test“ in the main menu to start the testing process. Here, you will receive an animated instruction video directly on your smartphone.

Before the actual test process starts, you will be guided through several steps for preparation (see **General Information**, page 67 and **Storage, Stability and Disposal**, page 69).

Fold open the paper stool catcher and pull it apart in the direction of the arrows. Carefully attach the paper stool catcher in the rear side of the toilet bowl using the **adhesive strips**.

Place the stool sample onto the paper stool catcher.

Attention: the stool sample should not come into contact with **toilet water or urine** before the sample is extracted.

To open the sample collection tube (TUBE), twist the yellow lid 90° and pull out the sampling stick.

Attention: Hold upright while opening to keep the buffer solution from leaking. Gently handle the tip of the TUBE to keep it from breaking.


Prick the stool sample in at least 3 different places, immediately one after the other, using the lower end of the yellow sampling stick.

Attention: All grooves of the sampling stick on both sides must be completely filled with stool.

Carefully remove stool from the blunt end of the sampling stick by wiping it on the paper stool catcher. The paper stool catcher containing the remaining stool can now be disposed of by flushing the toilet.

6

7

Insert the sampling stick with the stool sample into the TUBE filled with the buffering solution only **once** and **do not pull it back out**. Close the TUBE tightly by twisting the sampling stick 90°.


Shake the closed TUBE vigorously until the stool sample has been removed from the sampling stick and all grooves are free from stool.
Note: This causes the buffer solution in the TUBE to undergo only a slight color change.

Test Procedure

1

2

Remove the test cassette from the packaging and place it on a flat, dry and light surface.


Attention: Ensure that the test cassette lies flat at all times. Immediately carry out the next step.


Hold the TUBE with the top up, wrap it into soft paper (i.e. toilet paper) and then break the tip carefully.

Hold the TUBE approx. 1 cm above the sample application window. Place 3 hanging drops onto the testing cassette by pushing the TUBE in the middle at the height of the label.

Attention: the drops must freely fall from approximately 1 cm, to ensure the correct drop volume.

4**5**

Immediately after adding the sample, push „Start Timer“ in the SmarTest Calprotectin app. A 15-minute timer starts. A signal tone sounds after the timer expires.


You have **2 minutes** for evaluation. Align the stencil on the display with the test cassette — hold your smartphone parallel to the test cassette. **Attention:** Make sure that the lighting is even. Avoid shadows and direct sunlight. If the time is exceeded, the result becomes **invalid**. The measurement **must be repeated** with a new test kit.

6

The result is calculated by the app (using a 4-parameter fit) and is **immediately displayed on your smartphone**. The „questionnaire“ and „Explanation of results“ functions provide access to additional important medical information regarding your test results in the app. In some cases, you may be prompted to contact your physician. Depending on the settings, your results will be **transmitted automatically to the medical professionals involved in your care**. These settings can be configured under „mySmarTest“ on the main screen at any time (see **Releasing test results**, page 72).

Evaluation

The test result is **valid**, if:

- control line (C) appears.
- control line (C) appears and one or two test lines (T1 and/or T2).

The test result is **invalid**, if:

- **no** control line (C) appears.
- **no** control line (C) appears, but instead one or two test lines (T1 and/or T2) appears.

Test performance

Analytical performance

- The test can detect calprotectin concentrations between 31 and 4000 µg calprotectin per g stool.
- 22 different drugs (e.g. azathioprine and mesalazine) and food supplements (e.g. 25-OH vitamin D₃ and Ferro Sanol) were tested to determine whether the test reacts to them (cross-reaction). No relevant cross-reaction could be detected for any of the tested substances. A detailed list is available on request.
- A comparative measurement of 101 patient samples with the reference method ELISA (enzyme-linked immunosorbent assay) – a measurement method used in clinical laboratories – showed very good correlation between the two methods (89.5 %).
- Accuracy, repeatability and reproducibility tests were performed with different test batches, smartphones and in different environments. Variations were always smaller than 35 %.

Clinical performance

- A test result of ≥ 250 µg calprotectin per g stool is indicative of inflammatory events in the intestine. Healthy children under 4 years of age often have elevated calprotectin levels. In ulcerative colitis (UC) and Crohn's disease (CD), the measured calprotectin level provides information about the severity of the inflammatory event and can therefore indicate and even predict relapse.
- The measurement of calprotectin levels in stool can indicate active inflammatory disease in the intestine (e.g. UC or CD) with up to 81 % clinical accuracy and precision.
- Measurement of calprotectin in stool agrees very well (≥ 86 %) with various clinical scores (e.g. Mayo Score S-MES, Crohn's Disease Score).

Limitations of the test

- Although **SmarTest Calprotectin Home** reliably detects calprotectin, false results may occur in individual cases.
- For inconclusive results, we recommend contacting your medical professional to clarify the pathology and implement therapeutic measures.
- If results are questionable, it is recommended that a medical professional conducts additional clinical tests to finalize the diagnosis. As is the case for all diagnostic tests, the final diagnosis should not be based on a single result, but on a comprehensive clarification of the pathology by the medical professional responsible for treatment. The use of **SmarTest Calprotectin Home** is an additional tool for the optimum monitoring of the therapy for and progression of the disease and does not replace a medical professional's appointment.

References

1. Kawashima K, Ishihara S, Yuki T et al. Fecal calprotectin level correlated with both endoscopic severity and disease extent in ulcerative colitis. *BMC Gastroenterology* (2016) 16:47
2. Najari A., Hiraoka S, Takahashi S. et al. Simultaneous Measurements of Faecal Calprotectin and the Faecal Immunochemical Test in Quiescent Ulcerative Colitis Patients Can Stratify Risk of Relapse. *Journal of Crohn's and Colitis* (2018), 71–76
3. Amcoff K, Stridsberg M, Lampinen M et al. Clinical implications of assay specific differences in f-calprotectin when monitoring inflammatory bowel disease activity over time. *Scandinavian Journal of Gastroenterology* (2017), 52:3, 344-350
4. Matsura R., Watanabe O, Nakamura M. et al. Fecal calprotectin reflects endoscopic activity in patients with small-bowel Crohn's disease according to double-balloon endoscopy findings. *Nagoya J. Med. Sci.* 80 (2018). 257–266

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Literatür

1. Kawashima K, Ishihara S, Yuki T et al. Fecal calprotectin level correlated with both endoscopic severity and disease extent in ulcerative colitis. *BMC Gastroenterology* (2016) 16:47
2. Najari A, Hiraoka S, Takahashi S. et al. Simultaneous Measurements of Faecal Calprotectin and the Faecal Immunochemical Test in Quiescent Ulcerative Colitis Patients Can Stratify Risk of Relapse. *Journal of Crohn's and Colitis* (2018), 71–76
3. Amcoff K, Stridsberg M, Lampinen M et al. Clinical implications of assay specific differences in f-calprotectin when monitoring inflammatory bowel disease activity over time. *Scandinavian Journal of Gastroenterology* (2017), 52:3, 344-350
4. Matsuura R., Watanabe O, Nakamura M. et al. Fecal calprotectin reflects endoscopic activity in patients with small-bowel Crohn's disease according to double-balloon endoscopy findings. *Nagoya J. Med. Sci.* 80 (2018). 257–266

Immundiagnostik AG
Stubenwald-Allee 8a
64625 Bensheim, Germany

Distributie/Distribution/Vertrieb/
Distribution/Dağıtım:

Tramedico BV
Korte Muiderweg 2
1382 LR Weesp, Nederland

2021-11-23